



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RAB7 ANTIBODY

SKU: EB05479

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Ras-associated protein RAB7|PRO2706|RAB7, member RAS oncogene family|RAB7

Names

Peptide

Accession NP 004628 ID

Blocking EBP05479

Immunogen Peptide with sequence SFSSRPVPEILKT-C, from the N Terminus of the protein sequence according to NP 004628.

Peptide SFSSRPVPEILKT-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Refrigerated

Instructions

Predicted Human

Species

Human 7879

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western Blot: Preliminary experiments gave no signal but low background in Human Brain and Placenta lysates at Western **Blot**

up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with

other antibodies/lysates?