

GOAT ANTI-RAB7 ANTIBODY

SKU: EB05479



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Ras-associated protein RAB7 PRO2706 RAB7, member RAS oncogene family RAB7
Accession ID	NP_004628
Blocking Peptide	EBP05479
Immunogen	Peptide with sequence SFSSRPVPEILKT-C, from the N Terminus of the protein sequence according to NP_004628.
Peptide Sequence	SFSSRPVPEILKT-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	7879
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western Blot: Preliminary experiments gave no signal but low background in Human Brain and Placenta lysates at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?