



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RAB8A ANTIBODY

SKU: EB08311



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / ras-associated protein RAB8|mel transforming oncogene (derived from cell line NK14)- RAB8 homolog|mel transforming oncogene (derived from cell line NK14)|mel transforming oncogene (RAB8 homolog)|mel

transforming oncogene|RAB8|MEL|RAB8A, member RAS oncogene family|RAB8A Names

Accession

Alias

ID

NP_005361.2

Blocking

EBP08311 **Peptide**

Immunogen

Peptide with sequence C-NQGVKITPDQQKR, from the C Terminus of the protein sequence according to

NP 005361.2.

Peptide

C-NOGVKITPDOOKR Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Dog, Cow **Species**

Reactive

Human **Species**

Human

4218 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade

ELISA

Antibody detection limit dilution 1:32000. **Detection**

Limit

Western Approx 25kDa band observed in Human Brain (Cerebellum), Kidney and Heart lysates (calculated MW of

23.7kDa according to NP_005361.2). Recommended concentration: 0.1-0.3µg/ml. **Blot**

Application

Type

Pep-ELISA, WB

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa