

GOAT ANTI-RAB8A ANTIBODY

SKU: EB08311



250kDa
150kDa
100kDa
75kDa

50kDa

37kDa

25kDa
20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	ras-associated protein RAB8 mel transforming oncogene (derived from cell line NK14)- RAB8 homolog mel transforming oncogene (derived from cell line NK14) mel transforming oncogene (RAB8 homolog) mel transforming oncogene RAB8 MEL RAB8A, member RAS oncogene family RAB8A
Accession ID	NP_005361.2
Blocking Peptide	EBP08311
Immunogen	Peptide with sequence C-NQGVKITPDQQR, from the C Terminus of the protein sequence according to NP_005361.2.
Peptide Sequence	C-NQGVKITPDQQR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog, Cow
Reactive Species	Human
Human Gene ID	4218
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 25kDa band observed in Human Brain (Cerebellum), Kidney and Heart lysates (calculated MW of 23.7kDa according to NP_005361.2). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

