


GOAT ANTI-RABENOSYN 5 ANTIBODY

SKU: EB06053



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa

25kDa
20kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Rabenosyn-5 MGC126210 FLJ34993 zinc finger, FYVE domain containing 20 zinc finger, FYVE domain containing 20 FYVE-finger-containing Rab5 effector protein rabenosyn-5 Rabenosyn 5 ZFYVE20
Accession ID	NP_071735.2
Blocking Peptide	EBP06053
Immunogen	Peptide with sequence C-ELKHTLAKQKGGTD, from the C Terminus of the protein sequence according to NP_071735.2.
Peptide Sequence	C-ELKHTLAKQKGGTD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	64145
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 110 kDa band observed in A431, Fibroblast (epidermoid carcinoma) lysates (calculated MW of 88.8 kDa according to NP_071735). Band migration is consistent with literature (PMID 11062261, 15020713). Recommended concentration: 0.3 - 1.5 µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

