

GOAT ANTI-RABIF ANTIBODY

SKU: EB08547



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	mammalian suppressor of SEC4 guanine nucleotide exchange factor MSS4 Ras-specific guanine-releasing factor 3 RAB-interacting factor RASGRF3 RASGFR3 MSS4 RAB interacting factor RABIF
Accession ID	NP_002862.2
Blocking Peptide	EBP08547
Immunogen	Peptide with sequence C-ENVGFTKDVGNIK, from the internal region of the protein sequence according to NP_002862.2.
Peptide Sequence	C-ENVGFTKDVGNIK
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Cow
Reactive Species	Human
Human Gene ID	5877
Mouse Gene ID	98710
Rat Gene ID	304807
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx. 17kDa band observed in Human Prostate and Spleen lysates (calculated MW of 13.8kDa according to NP_002862.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

