



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RABIF ANTIBODY

SKU: EB08547



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Unit Size

100 µg

Synonym / Alias

mammalian suppressor of SEC4|guanine nucleotide exchange factor MSS4|Ras-specific guanine-releasing

factor 3|RAB-interacting factor|RASGRF3|RASGFR3|MSS4|RAB interacting factor|RABIF

Names Accession

NP_002862.2

Blocking

EBP08547

Peptide Immunogen

Peptide with sequence C-ENVGFTKDVGNIK, from the internal region of the protein sequence according to

NP 002862.2.

Peptide Sequence

C-ENVGFTKDVGNIK

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat, Cow

Reactive

Species

Human

Human Gene ID

5877

Mouse

98710

Gene ID

Rat Gene ID 304807

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade **ELISA**

Antibody detection limit dilution 1:2000. Detection

Limit

Western Approx. 17kDa band observed in Human Prostate and Spleen lysates (calculated MW of 13.8kDa according to

Blot NP_002862.2). Recommended concentration: 1-3µg/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa