

GOAT ANTI-RAC1 ANTIBODY

SKU: EB05306



SPECIFICATIONS

| | |
|------------------------------|---|
| Formulation | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. |
| Unit Size | 100 µg |
| Storage Instructions | Aliquot and store at -20°C. Minimize freezing and thawing. |
| Synonym / Alias Names | rho family, small GTP binding protein Rac1 TC-25 ras-related C3 botulinum toxin substrate 1 (rho family, small GTP binding protein Rac1) RAC1 |
| Accession ID | NP_008839 |
| Immunogen | Peptide with sequence PPPVKRKRKCLLL, from the C Terminus of the protein sequence according to NP_008839. |
| Product Comments | This antibody is expected to recognise both human isoforms of this protein. |
| Peptide Sequence | PPPVKRKRKCLLL |
| Purification Method | Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. |
| Shipping Instructions | Refrigerated |
| Predicted Species | Human, Mouse, Dog, Cow, Chicken |
| Reactive Species | Human |
| Human Gene ID | 5879 |
| Product Grade | https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png |
| Western Blot | Western blot: Approx 30kDa band seen in Human Brain lysate. Recommended for use at 2-4µg/ml. |