


GOAT ANTI-RANBP7 / IMPORTIN 7 ANTIBODY

SKU: EB05553



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	RAN-binding protein 7 RAN binding protein 7 MGC138673 FLJ14581 RANBP7 importin 7 IPO7
Usage Summary	Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration: 10ug/ml.
Accession ID	NP_006382.1
Blocking Peptide	EBP05553
Immunogen	Peptide with sequence C-SSFNFGGPAPGMN, from the C Terminus of the protein sequence according to NP_006382.1.
Peptide Sequence	C-SSFNFGGPAPGMN
Shipping Instructions	Refrigerated

Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	10527
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA, FC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

