


# GOAT ANTI-RAP80 ANTIBODY

**SKU:** EB09551



## SPECIFICATIONS

**Formulation** Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

**Storage Instructions** Aliquot and store at -20°C. Minimize freezing and thawing.

**Synonym / Alias Names** retinoid x receptor interacting protein|receptor associated protein 80|OTTHUMP00000161441|X2HRIP110|RAP80|ubiquitin interaction motif containing 1|UIMC1

**Accession ID** NP\_057374.3

**Blocking Peptide** EBP09551

**Immunogen** Peptide with sequence C-TRRQKEAKTKSDSG, from the internal region of the protein sequence according to NP\_057374.3.

**Peptide Sequence** C-TRRQKEAKTKSDSG

**Purification Method** Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

**Shipping Instructions** Refrigerated

<b>Predicted Species</b>	Human, Mouse, Rat, Dog, Pig
<b>Reactive Species</b>	Human
<b>Human Gene ID</b>	51720
<b>Mouse Gene ID</b>	20184
<b>Rat Gene ID</b>	290997
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>
<b>IHC Results</b>	Paraffin embedded Human Testis. Recommended concentration: 3.75µg/ml.
<b>ELISA</b>	
<b>Detection Limit</b>	Antibody detection limit dilution 1:32000.
<b>Application Type</b>	Pep-ELISA, IHC

## DOCUMENTS

- [Data Sheet](#)

## GALLERY IMAGES

