


Telephone: (650) 697-3600

## **GOAT ANTI-RAP80 ANTIBODY**

**SKU:** EB09551



## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

**Alias** 

retinoid x receptor interacting protein|receptor associated protein

80|OTTHUMP00000161441|X2HRIP110|RAP80|ubiquitin interaction motif containing 1|UIMC1 **Names** 

**Accession** 

ID

NP\_057374.3

**Blocking** 

EBP09551

**Peptide** 

Peptide with sequence C-TRRQKEAKTKSDSG, from the internal region of the protein sequence according to

NP\_057374.3.

**Peptide** 

Sequence

Immunogen

C-TRRQKEAKTKSDSG

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

**Predicted** 

Human, Mouse, Rat, Dog, Pig

**Species** Reactive

Human

**Species** Human

51720

Gene ID Mouse

**Gene ID** 

20184

**Rat Gene ID** 290997

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

IHC Results Paraffin embedded Human Testis. Recommended concentration: 3.75µg/ml.

**ELISA** 

Detection

Antibody detection limit dilution 1:32000.

Limit

**Type** 

**Application** Pep-ELISA, IHC

## **GALLERY IMAGES**

