

# EBP06673 - IMMUNIZING PEPTIDE

SKU: EBP06673



## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Shipped at ambient temperature, store at -20°C.
<b>Accession ID</b>	NP_000956.2; NP_057236.1
<b>Immunogen</b>	Peptide with sequence C-VENSGVSQSPLVQ, from the C Terminus of the protein sequence according to NP_000956.2; NP_057236.1.
<b>Product Comments</b>	The peptide was used in the production of EB06673 - Goat Anti-RARB Antibody. It is currently untested in blocking applications.
<b>Peptide Sequence</b>	C-VENSGVSQSPLVQ
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>
<b>Western Blot</b>	Preliminary experiments in Human Brain, Bone Marrow and Liver lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

## DOCUMENTS

- [Data Sheet](#)