

GOAT ANTI-RARS ANTIBODY

SKU: EB09780



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names RARS|MGC8641|DALRD1|ArgRS|arginyl-tRNA synthetase|arginine tRNA ligase 1, cytoplasmic

Accession ID NP_002878.2

Blocking Peptide EBP09780

Immunogen Peptide with sequence C-SLQAERNKPTKN, from the internal region of the protein sequence according to NP_002878.2.

Peptide Sequence C-SLQAERNKPTKN

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat, Cow

Reactive Species Human

Human Gene ID 5917

Mouse Gene ID 104458

Rat Gene ID 287191

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Limit Antibody detection limit dilution 1:128000.

Western Blot Approx 75kDa band observed in lysates of HeLa and HepG2 (calculated MW of 75.4kDa according to NP_002878.2). Recommended concentration: 1-3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

