

GOAT ANTI-RASAL ANTIBODY

SKU: EB06176



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa
10kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	RAS protein activator like 1 RAS-GTPase-activating protein-like GAP1 like protein RAS protein activator like 1 (GAP1 like) RASAL RASAL1
Accession ID	NP_004649.1
Blocking Peptide	EBP06176
Immunogen	Peptide with sequence C-QERGKAALGPLGP, from the C Terminus of the protein sequence according to NP_004649.1.
Peptide Sequence	C-QERGKAALGPLGP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	8437
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 90kDa band observed in Human Thyroid gland lysates (calculated MW of 90.0kDa according to NP_004649.1). Recommended concentration: 0.3-1µg/ml. An additional band of unknown identity was also consistently observed at 21kDa. This band was successfully blocked by incubation with the immunising peptide.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

