



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RASAL ANTIBODY

SKU: EB06176





Telephone: (650) 697-3600

250kDa 150kDa 100kDa

e erest

75kDa

50kDa


37kDa

25kDa

20kDa

15kDa

10kDa





SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

RAS protein activator like 1|RAS-GTPase-activating protein-like|GAP1 like protein|RAS protein activator like 1

(GAP1 like)|RASAL|RASAL1

Names **Accession**

ID

NP_004649.1

Blocking

Peptide

EBP06176

Peptide with sequence C-QERGKAALGPLGP, from the C Terminus of the protein sequence according to **Immunogen**

NP 004649.1.

Peptide

Sequence

C-QERGKAALGPLGP

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human

Human

Reactive

Species

Human

8437 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Western

Antibody detection limit dilution 1:32000.

Limit

Approx 90kDa band observed in Human Thyroid gland lysates (calculated MW of 90.0kDa according to NP_004649.1). Recommended concentration: 0.3-1μg/ml. An additional band of unknown identity was also consistently observed at 21kDa. This band was successfully blocked by incubation with the immunising

peptide.

Application

Type

Blot


Pep-ELISA, WB

GALLERY IMAGES





Telephone: (650) 697-3600



e erest

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa