

GOAT ANTI-RASSF1 ANTIBODY

SKU: EB11904



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	tumor suppressor protein RDA32 ras association domain-containing protein 1 pancreas-specific ras association domain family 1 protein cardiac-specific ras association domain family 1
Names	protein WUGSC:H_LUCA12.5 OTTHUMP00000195507 OTTHUMP00000195505 OTTHUMP00000195504 OTTHUMP00000195503 REH3P21 RDA32 RASSF1A NORE2A 123F2 Ras association (RalGDS/AF-6) domain family member 1 RASSF1
Accession ID	NP_009113.3; NP_733830.1; NP_733831.1; NP_733832.1
Blocking Peptide	EBP11904
Immunogen	Peptide with sequence C-RAREVIEALLRKFLVDDPRK, from the C Terminus of the protein sequence according to NP_009113.3; NP_733830.1; NP_733831.1; NP_733832.1.
Product Comments	This antibody is expected to recognize all reported isoforms (NP_009113.3; NP_733830.1; NP_733831.1; NP_733832.1). Reported variants represent identical protein: NP_733830.1, NP_001193886.1
Peptide Sequence	C-RAREVIEALLRKFLVDDPRK
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog, Pig, Cow
Reactive Species	Human
Human Gene ID	11186
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 40kDa band observed in Human Heart lysates (calculated MW of 38.8kDa according to NP_009113.3). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

