

GOAT ANTI-RASSF1 ANTIBODY

SKU: EB11421



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	RASSF1 WUGSC:H_LUCA12.5 tumor suppressor protein RDA32 REH3P21 RDA32 RASSF1A ras association domain-containing protein 1 Ras association (RalGDS/AF-6) domain family member 1 pancreas-specific ras association domain family 1 protein OTTHUMP00000195507 OTTHUMP00000195505
Names	OTTHUMP00000195504 OTTHUMP00000195503 NORE2A cardiac-specific ras association domain family 1 protein 123F2
Accession ID	NP_009113.3; NP_733830.1; NP_733831.1; NP_733832.1
Blocking Peptide	EBP11421
Immunogen	Peptide with sequence QKYSYCRQKIQEALH, from the C Terminus of the protein sequence according to NP_009113.3; NP_733830.1; NP_733831.1; NP_733832.1.
Product Comments	This antibody is expected to recognize all reported isoforms (NP_009113.3; NP_733830.1; NP_733831.1; NP_733832.1). Reported variants represent identical protein: NP_733830.1, NP_001193886.1
Peptide Sequence	QKYSYCRQKIQEALH
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog, Pig, Cow
Reactive Species	Human
Human Gene ID	11186
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 38kDa band observed in lysates of cell lines A549 and MRC5 PD19 (calculated MW of 38.8kDa according to NP_009113.3). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

