
Telephone: (650) 697-3600





GOAT ANTI-RASSF2 ANTIBODY

SKU: EB08511



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias

rasfadin|Ras association domain family 2|DKFZp78101747N: KIAA0168|Ras association (RalGDS/AF-6) domain

family 2|RASSF2 **Names**

Accession

ID

NP_055552.1; NP_739580.1

Blocking Peptide

EBP08511

Immunogen

Peptide with sequence PCGQDKYISKNE, from the N Terminus of the protein sequence according to

NP_055552.1; NP_739580.1.

Product Comments

Reported variants represent identical protein: NP 055552.1; NP 739580.1.

Peptide Sequence

PCGQDKYISKNE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Shipping

Refrigerated Instructions

Predicted

Human, Cow **Species**

Reactive

Human **Species**

e erest

Human

9770

Gene ID **Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results Paraffin embedded Human Spleen, Thymus and Liver. Recommended concentration: 3.75µg/ml.

ELISA

Detection Antibody detection limit dilution 1:16000.

Limit

Type

Application Pep-ELISA, IHC

SELECTED REFERENCES

[{"pmid": 26284587, "intro": "This antibody has been successfully used in IHC and IF in Human:", "title":

"Differential involvement of RASSF2 hypermethylation in breast cancer subtypes and their prognosis.", "author": "Perez-Janices N, Blanco-Luquin I, Torrea N, Liechtenstein T, Escors D, Cordoba A, Vicente-Garcia F, Jauregui I, De La Cruz S, Illarramendi JJ, Coca V, Berdasco M, Kochan G, Ibañez B, Lera JM, Guerrero-Setas D.", "journal": "Oncotarget. 2015 Sep 15;6(27):23944-58."}]

GALLERY IMAGES

