




GOAT ANTI-RASSF2 ANTIBODY

SKU: EB08511



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	rasfadin Ras association domain family 2 DKFZp781O1747N: KIAA0168 Ras association (RalGDS/AF-6) domain family 2 RASSF2
Accession ID	NP_055552.1; NP_739580.1
Blocking Peptide	EBP08511
Immunogen	Peptide with sequence PCGQDKYISKNE, from the N Terminus of the protein sequence according to NP_055552.1; NP_739580.1.
Product Comments	Reported variants represent identical protein: NP_055552.1; NP_739580.1.
Peptide Sequence	PCGQDKYISKNE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.



Shipping Instructions	Refrigerated
Predicted Species	Human, Cow
Reactive Species	Human
Human Gene ID	9770
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Spleen, Thymus and Liver. Recommended concentration: 3.75µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA, IHC

SELECTED REFERENCES

[{"pmid": 26284587, "intro": "**This antibody has been successfully used in IHC and IF in Human:**", "title": "Differential involvement of RASSF2 hypermethylation in breast cancer subtypes and their prognosis.", "author": "Perez-Janices N, Blanco-Luquin I, Torrea N, Liechtenstein T, Escors D, Cordoba A, Vicente-Garcia F, Jauregui I, De La Cruz S, Illarramendi JJ, Coca V, Berdasco M, Kochan G, Ibañez B, Lera JM, Guerrero-Setas D.", "journal": "Oncotarget. 2015 Sep 15;6(27):23944-58."}]

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

