

GOAT ANTI-RASSF3 ANTIBODY

SKU: EB08391



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names RASSF5|MGC119197|MGC119195|MGC119194|Ras association (RalGDS/AF-6) domain family 3|RASSF3

Accession ID NP_835463.1

Blocking Peptide EBP08391

Immunogen Peptide with sequence CKPPQTSPNSGK, from the C Terminus of the protein sequence according to NP_835463.1.

Peptide Sequence CKPPQTSPNSGK

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human

Reactive Species Human

Human Gene ID 283349

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Limit Antibody detection limit dilution 1:2000.

Western Blot Approx 26kDa band observed in Human Heart lysates (calculated MW of 27.6kDa according to NP_835463.1). Recommended concentration: 0.5-1.5µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

