



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RASSF3 ANTIBODY

SKU: EB08391

Telephone: (650) 697-3600






Email: customerservice@vectorlabs.com


250kDa 150kDa 100kDa 75kDa 50kDa

37kDa

25kDa 20kDa

15kDa







SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

RASSF5|MGC119197|MGC119195|MGC119194|Ras association (RaIGDS/AF-6) domain family 3|RASSF3

Names

ID

Accession

NP_835463.1

Blocking

Peptide

EBP08391

Immunogen

Peptide with sequence CKPPQTSPNSGK, from the C Terminus of the protein sequence according to

NP 835463.1.

Peptide

Sequence

CKPPQTSPNSGK

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

Reactive

Species

Human

283349

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade

ELISA

Detection

Antibody detection limit dilution 1:2000.

Limit

Western Approx 26kDa band observed in Human Heart lysates (calculated MW of 27.6kDa according to NP_835463.1).

Blot Recommended concentration: 0.5-1.5µg/ml.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES





Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa