

GOAT ANTI-RASSF4 ANTIBODY

SKU: EB08392



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	tumor suppressor RASSF4 Ras association domain family 4 MGC44914 AD037 RP11-285G1.4 Ras association (RalGDS/AF-6) domain family 4 RASSF4
Accession ID	NP_114412.2
Blocking Peptide	EBP08392
Immunogen	Peptide with sequence HESGERTKLKDC, from the internal region of the protein sequence according to NP_114412.2.
Peptide Sequence	HESGERTKLKDC
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Cow
Human Gene ID	83937
Mouse Gene ID	213391
Rat Gene ID	362423
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)