Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RASSF4 ANTIBODY

SKU: EB08392

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

tumor suppressor RASSF4|Ras association domain family 4|MGC44914|AD037|RP11-285G1.4|Ras association **Alias** (RaIGDS/AF-6) domain family 4|RASSF4

Names

Accession

NP 114412.2 ID

Blocking

EBP08392 **Peptide**

Peptide with sequence HESGERTKLKDC, from the internal region of the protein sequence according to Immunogen

NP 114412.2.

Peptide

HESGERTKLKDC Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Cow **Species**

Human

Gene ID

83937

Mouse **Gene ID**

213391

Rat Gene ID 362423

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:4000.

Limit

Application Pep-ELISA

Type