

GOAT ANTI-RASSF6 ANTIBODY

SKU: EB08394



250kDa
150kDa
100kDa

75kDa
50kDa

37kDa

25kDa
20kDa

15kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names RASSF6|Ras association (RalGDS/AF-6) domain family 6|DKFZp686K23225|putative RAS binding protein

Accession ID NP_803876.1; NP_958834.1

Blocking Peptide EBP08394

Immunogen Peptide with sequence C-SNTLKPHAKDEPDS, from the internal region of the protein sequence according to NP_803876.1; NP_958834.1.

Product Comments This antibody is expected to recognise both reported isoforms (NP_803876.1; NP_958834.1).

Peptide Sequence C-SNTLKPHAKDEPDS

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human

Reactive Species Human

Human Gene ID 166824

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:16000.

Western Blot Approx 45kDa band observed in lysates of cell line HepG2 (calculated MW of 43.4kDa according to NP_958834.1). Recommended concentration: 1-3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

