



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RASSF6 ANTIBODY

SKU: EB08394



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa 50kDa

37kDa

25kDa 20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / **Alias** RASSF6|Ras association (RalGDS/AF-6) domain family 6|DKFZp686K23225|putative RAS binding protein

Names

Accession

NP_803876.1; NP_958834.1

Blocking

EBP08394 **Peptide**

Peptide with sequence C-SNTLKPHAKDEPDS, from the internal region of the protein sequence according to **Immunogen**

NP_803876.1; NP_958834.1.

Product Comments

This antibody is expected to recognise both reported isoforms (NP 803876.1; NP 958834.1).

Peptide

Sequence

C-SNTLKPHAKDEPDS

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive **Species**

Human

Human

166824 Gene ID

Product

Grade **ELISA**

Detection Antibody detection limit dilution 1:16000.

Limit

Western Approx 45kDa band observed in lysates of cell line HepG2 (calculated MW of 43.4kDa according to

NP_958834.1). Recommended concentration: 1-3µg/ml. Blot

Application

Pep-ELISA, WB Type

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa 50kDa

37kDa

25kDa 20kDa

15kDa