

GOAT ANTI-RASSF7 / HRAS1 ANTIBODY

SKU: EB08383



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	HRAS1-related cluster-1 MGC126070 MGC126069 HRC1 C11orf13 Ras association (RalGDS/AF-6) domain family 7 HRAS1 RASSF7
Accession ID	NP_003466.1
Blocking Peptide	EBP08383
Immunogen	Peptide with sequence CTDLRGLELRVQRN, from the internal region of the protein sequence according to NP_003466.1.
Peptide Sequence	CTDLRGLELRVQRN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	8045
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 37kDa band observed in lysates of cell line Jurkat (calculated MW of 39.9kDa according to NP_003466.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

