

GOAT ANTI-RASSF8 ANTIBODY

SKU: EB08395



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	RASSF8 Ras association (RalGDS/AF-6) domain family 8 C12orf2 HoJ-1 carcinoma associated HOJ-1
Accession ID	NP_009142.2
Blocking Peptide	EBP08395
Immunogen	Peptide with sequence C-RDTERHLAPHEN, from the internal region of the protein sequence according to NP_009142.2.
Peptide Sequence	C-RDTERHLAPHEN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Cow
Reactive Species	Human
Human Gene ID	11228
Mouse Gene ID	71323
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 40-45kDa band observed in Human Testis lysates (calculated MW of 45.4kDa according to NP_009142.2). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

