



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RASSF8 ANTIBODY

SKU: EB08395



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa


75kDa

50kDa

37kDa

25kDa 20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Instructions Aliquot and store at -20°C. Minimize freezing and thawing. Storage

Synonym /

Alias Names

Accession NP_009142.2

Blocking EBP08395 **Peptide**

Peptide with sequence C-RDTERHLAPHEN, from the internal region of the protein sequence according to **Immunogen**

RASSF8|Ras association (RalGDS/AF-6) domain family 8|C12orf2|HoJ-1|carcinoma associated HOJ-1

NP 009142.2.

Peptide C-RDTERHLAPHEN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Refrigerated Instructions

Predicted

Human, Mouse, Cow **Species**

Reactive

Human **Species**

Human

11228 Gene ID

Mouse 71323

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:128000.

Limit

Western Approx 40-45kDa band observed in Human Testis lysates (calculated MW of 45.4kDa according to

Blot NP_009142.2). Recommended concentration: 0.3-1µg/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES





Telephone: (650) 697-3600

e e e rest

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa 20kDa

15kDa