

GOAT ANTI-RBP1-LIKE PROTEIN ANTIBODY

SKU: EB05165



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	breast carcinoma-associated antigen BCAA retinoblastoma-binding protein 1-like 1 breast cancer-associated antigen BRCA1 Rb-binding protein homolog RBBP1L1 BRCA1 RBP1-like protein BCAA
Accession ID	NP_057458
Immunogen	Peptide with sequence SSASQNGMSVECR, from the C Terminus of the protein sequence according to NP_057458.
Product Comments	This antibody will recognise both isoforms.
Peptide Sequence	SSASQNGMSVECR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	51742
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: Two consistent bands observed in HeLa at approx 45kDa and approx 65kDa with 12 hour incubation of primary.