

GOAT ANTI-RE1-SILENCING TRANSCRIPTION FACTOR ANTIBODY

SKU: EB12701



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	repressor binding to the X2 box neuron restrictive silencer factor neural-restrictive silencer factor RE1-silencing transcription factor variant E1c/E2g/E3/E4 RE1-silencing transcription factor variant E1c/E2/E3/E5 RE1-silencing transcription factor variant E1b/E2c/E2j/E3/E4 RE1-silencing transcription factor variant E1b/E2/E3/N3c/E4 RE1-silencing transcription factor variant E1b/E2/E3/N3b/E4i RE1-silencing transcription factor variant E1b/E2/E3/E5 RE1-silencing transcription factor variant E1a/E2k/E2i/E3/E4j RE1-silencing transcription factor variant E1a/E2/E4 RE1-silencing transcription factor variant E1a/E2/E3/N3c/E4 RE1-silencing transcription factor variant E1a/E2/E3/N3a/E4i RE1-silencing transcription factor variant E1a/E2/E3/E5 RE1-silencing transcription factor variant E1a/E2/E3/E4c XBR NRSF RE1-silencing transcription factor REST
Accession ID	NP_005603.3
Blocking Peptide	EBP12701
Immunogen	Peptide with sequence C-PETPGAEDKGKSSK, from the internal region of the protein sequence according to NP_005603.3.
Product Comments	Reported variants represent identical protein: NP_001180437.1, NP_005603.3
Peptide Sequence	C-PETPGAEDKGKSSK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	5978
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 140kDa band observed in nuclear lysates of cell lines HEK293 and HeLa (calculated MW of 122kDa according to NP_005603.3). This molecular weight is routinely observed by other sources. Recommended concentration: 0.3-1µg/ml.

Application
Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES



