

GOAT ANTI-RENT2 / UPF2 ANTIBODY

SKU: EB07075



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	OTTHUMP00000019091 smg-3 MGC138835 MGC138834 yeast Upf2p homolog regulator of nonsense transcripts 2 UPF2 regulator of nonsense transcripts homolog OTTHUMP00000045025 RENT2 KIAA1408 HUPF2 DKFZP434D222 HGNC:17854 UPF2 regulator of nonsense transcripts homolog (yeast) UPF2
Accession ID	NP_542166.1; NP_056357.1
Blocking Peptide	EBP07075
Immunogen	Peptide with sequence C-PRYQHPKGAPNAD, from the C Terminus of the protein sequence according to NP_542166.1; NP_056357.1.
Product Comments	This antibody is expected to recognise both reported isoforms (NP_542166.1; NP_056357.1).
Peptide Sequence	C-PRYQHPKGAPNAD
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	26019
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Preliminary experiments in Human HeLa lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)