Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RENT2 / UPF2 ANTIBODY

SKU: EB07075

SPECIFICATIONS

Unit Size 100 µg

Synonym /

OTTHUMP00000019091|smg-3|MGC138835|MGC138834|yeast Upf2p homolog|regulator of nonsense transcripts

Alias

2|UPF2 regulator of nonsense transcripts

Names

homolog|OTTHUMP00000045025|RENT2|KIAA1408|HUPF2|DKFZP434D222|HGNC:17854|UPF2 regulator of

nonsense transcripts homolog (yeast)|UPF2

Accession

ID

NP_542166.1; NP_056357.1

Blocking

Peptide

EBP07075

Peptide with sequence C-PRYQHPKGAPNAD, from the C Terminus of the protein sequence according to Immunogen NP_542166.1; NP_056357.1.

Product

Comments

This antibody is expected to recognise both reported isoforms (NP_542166.1; NP_056357.1).

Peptide

C-PRYQHPKGAPNAD Sequence

Refrigerated Instructions **Shipping**

Predicted

Human, Mouse, Rat, Dog, Cow

Species Human

26019 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection

Antibody detection limit dilution 1:4000.

Limit Western

Preliminary experiments in Human HeLa lysates gave no specific signal but low background (at antibody

concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been

reported with other antibodies/lysates?

Application

Type

Blot

Pep-ELISA