

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RERG ANTIBODY

SKU: EB05280

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias MGC15754|RAS-like, estrogen-regulated, growth-inhibitor|RERG

Names

Blocking EBP05280 **Peptide**

Immunogen Peptide with sequence C-VKQAINKMLTKISS, from the C Terminus of the protein sequence according to .

Peptide C-VKQAINKMLTKISS

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Refrigerated

Instructions

Predicted Human, Mouse

Human

Species

Sequence

Reactive

Species

Human Gene ID

85004

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Western

Western blot: Approx 23kDa band seen in 293 lysates. Recommended for use at 1-3µg/ml. Blot