

GOAT ANTI-RERG ANTIBODY

SKU: EB05280

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	MGC15754 RAS-like, estrogen-regulated, growth-inhibitor RERG
Names	
Blocking Peptide	EBP05280
Immunogen	Peptide with sequence C-VKQAINKMLTKISS, from the C Terminus of the protein sequence according to .
Peptide Sequence	C-VKQAINKMLTKISS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	85004
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 23kDa band seen in 293 lysates. Recommended for use at 1-3µg/ml.