

# GOAT ANTI-RGS16 ANTIBODY

**SKU:** EB05423



## SPECIFICATIONS

|                              |                                                                                                                                                                                                                                         |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>             | 100 µg                                                                                                                                                                                                                                  |
| <b>Synonym /</b>             |                                                                                                                                                                                                                                         |
| <b>Alias</b>                 | Regulator of G protein signaling-16 A28-RGS14 RGS-r regulator of G-protein signalling 16 RGS16                                                                                                                                          |
| <b>Names</b>                 |                                                                                                                                                                                                                                         |
| <b>Accession ID</b>          | NP_002919                                                                                                                                                                                                                               |
| <b>Blocking Peptide</b>      | EBP05423                                                                                                                                                                                                                                |
| <b>Immunogen</b>             | Peptide with sequence RTLAAPPTCLER, from the N Terminus of the protein sequence according to NP_002919.                                                                                                                                 |
| <b>Peptide Sequence</b>      | RTLAAPPTCLER                                                                                                                                                                                                                            |
| <b>Shipping Instructions</b> | Refrigerated                                                                                                                                                                                                                            |
| <b>Predicted Species</b>     | Human                                                                                                                                                                                                                                   |
| <b>Reactive Species</b>      | Human                                                                                                                                                                                                                                   |
| <b>Human Gene ID</b>         | 6004                                                                                                                                                                                                                                    |
| <b>Mouse Gene ID</b>         | 19734                                                                                                                                                                                                                                   |
| <b>Product Grade</b>         | <a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a> |
| <b>Western Blot</b>          | Western blot: Consistent strong band observed in Jurkat, Human Brain and Human Lung lysates at approx 40kDa.                                                                                                                            |