

GOAT ANTI-RGS17 / RGSZ2 ANTIBODY

SKU: EB06370



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	regulator of G-protein signalling 17 OTTHUMP00000017460 hRGS17 regulator of G-protein signalling 17 RGS-17 RGSZ2 RGS17
Accession ID	NP_036551.3
Blocking Peptide	EBP06370
Immunogen	Peptide with sequence C-SFVESTAGSSSES, from the C Terminus of the protein sequence according to NP_036551.3.
Peptide Sequence	C-SFVESTAGSSSES
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	26575
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)