

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RGS17 / RGSZ2 ANTIBODY

SKU: EB06370

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

regulator of G-protein signalling 17|OTTHUMP00000017460|hRGS17|regulator of G-protein signalling **Alias**

17|RGS-17|RGSZ2|RGS17 **Names**

Accession

NP 036551.3 ID

Blocking

EBP06370 **Peptide**

Peptide with sequence C-SFVESTAGSSSES, from the C Terminus of the protein sequence according to Immunogen

NP 036551.3.

Peptide

C-SFVESTAGSSSES Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human **Gene ID**

26575

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:16000.

Limit

Application Pep-ELISA **Type**