


Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-RHBG (MOUSE) ANTIBODY**

**SKU:** EB10280







250kDa 150kDa 100kDa 75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / RHBG|SLC42A2|Rhesus blood group, B glycoprotein|Rh type B glycoprotein|Rh family, B glycoprotein Alias

(gene/pseudogene) Names

**Accession** 

NP\_067350.2 ID

**Blocking** 

EBP10280 **Peptide** 

Peptide with sequence C-PLRGGESDTR, from the C Terminus of the protein sequence according to **Immunogen** 

NP 067350.2.

**Peptide** 

C-PLRGGESDTR Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Mouse **Species** 

Reactive

Mouse **Species** 

Mouse

58176 Gene ID

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:1000.

Limit

**Blot** 

Western Approx 55kDa band observed in Mouse Kidney lysates and in lysates of cell lines NIH3T3 and NSO (calculated

MW of 49.5kDa according to NP\_067350.2). Recommended concentration: 1-3µg/ml.

**Application** 

Pep-ELISA, WB Type

## **GALLERY IMAGES**



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: <u>(650)</u> 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa