


GOAT ANTI-RIG / DKK3 ANTIBODY

SKU: EB05940



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	DKK3 regulated in glioma RIG
Accession ID	AAC35996.1
Blocking Peptide	EBP05940
Immunogen	Peptide with sequence TLSVVGPHPKPCIL, from the C Terminus of the protein sequence according to AAC35996.1.
Product Comments	NCBI has removed RIG (GeneID 10530) and has it replaced by DKK3 (GeneID 27122). The AAC35996 protein is no longer represented by an NP RefSeq code.
Peptide Sequence	TLSVVGPHPKPCIL
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	10530
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Uterus shows membranous staining in the endometrial glands. Recommended concentration, 2-6µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

