




Telephone: (650) 697-3600



GOAT ANTI-RIG / DKK3 ANTIBODY

SKU: EB05940



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias

DKK3|regulated in glioma|RIG

Names

ID

Accession

AAC35996.1

Blocking

EBP05940

Peptide

Peptide with sequence TLSVVGPHPKPCIL, from the C Terminus of the protein sequence according to

Immunogen

Product

NCBI has removed RIG (GeneID 10530) and has it replaced by DKK3 (GeneID 27122). The AAC35996 protein is

Comments no longer represented by an NP RefSeq code.

Peptide Sequence

TLSVVGPHPKPCIL

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Shipping Instructions

Refrigerated

Predicted

Human **Species**

Reactive

Human **Species**

Human

10530

Gene ID Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$ Grade

IHC Results

In paraffin embedded Human Uterus shows membranous staining in the endometrial glands. Recommended

concentration, 2-6µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:16000.

Limit

Application

Type

Pep-ELISA, IHC

GALLERY IMAGES

