



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600



GOAT ANTI-RIN1 ANTIBODY

SKU: EB05219

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias ras inhibitor 1|Ras and Rab interactor 1|RIN1

Names

Usage Please note that the stronger 100 kDa band corresponds to the reference sequence NP 004283, whereas the

Summary fainter 35 kDa band may correspond to a novel isoform of this protein.

Accession

NP 004283

Peptide with sequence C-PGEPEAEGSRAAEE, from the C Terminus of the protein sequence according to **Immunogen**

NP 004283.

Peptide

C-PGEPEAEGSRAAEE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human

Species Human

9610

Gene ID **Product**

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

Western blot: Approx 100kDa band and 35kDa bands seen in HeLa cell lysates. Recommended for use at Western

Blot $1-3\mu g/ml$.