

GOAT ANTI-RIN1 ANTIBODY

SKU: EB05219



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	ras inhibitor 1 Ras and Rab interactor 1 RIN1
Names	
Usage Summary	Please note that the stronger 100 kDa band corresponds to the reference sequence NP_004283, whereas the fainter 35 kDa band may correspond to a novel isoform of this protein.
Accession ID	NP_004283
Immunogen	Peptide with sequence C-PGEPEAEGRSRAAEE, from the C Terminus of the protein sequence according to NP_004283.
Peptide Sequence	C-PGEPEAEGRSRAAEE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	9610
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 100kDa band and 35kDa bands seen in HeLa cell lysates. Recommended for use at 1-3µg/ml.