


GOAT ANTI-RIP140 ANTIBODY

SKU: EB09546



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / receptor interacting protein 140|nuclear factor

Alias RIP140|OTTHUMP00000083037|OTTHUMP00000083036|RIP140|FLJ77253|nuclear receptor interacting protein

Names 1|NRIP1

Accession ID NP_003480.2

Blocking Peptide EBP09546

Immunogen Peptide with sequence C-EKDSPRLTKTNP, from the internal region of the protein sequence according to NP_003480.2.

Peptide Sequence C-EKDSPRLTKTNP

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat, Dog, Cow

Reactive Species Human

Human Gene ID 8204

Mouse Gene ID 268903

Rat Gene ID 304157

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:32000.

Western Blot Approx 140kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 127kDa according to NP_003480.2). Recommended concentration: 0.5-2µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

