Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-RIPK2 ANTIBODY**

**SKU:** EB05390

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

**Alias** receptor interacting protein 2|CARDIAK|CARD3|RIP2|RICK|receptor-interacting serine-threonine kinase 2|RIPK2

**Names** 

**Accession** NP\_003812 ID

**Blocking** 

EBP05390 **Peptide** 

Peptide with sequence C-SRSPSLNLLQNKSM, from the C Terminus of the protein sequence according to **Immunogen** 

NP 003812.

**Peptide** C-SRSPSLNLLQNKSM

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** Refrigerated

Instructions

**Predicted** Human

**Species** 

Reactive Human **Species** 

Human

8767 Gene ID

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade

Western

Western blot: Approx 60-65kDa band seen in Daudi lysate. Recommended for use at 1-3μg/ml. **Blot**