

GOAT ANTI-RIPK2 ANTIBODY

SKU: EB05390



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	receptor interacting protein 2 CARDIAK CARD3 RIP2 RICK receptor-interacting serine-threonine kinase 2 RIPK2
Accession ID	NP_003812
Blocking Peptide	EBP05390
Immunogen	Peptide with sequence C-SRSPSLNLLQNKSM, from the C Terminus of the protein sequence according to NP_003812.
Peptide Sequence	C-SRSPSLNLLQNKSM
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	8767
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 60-65kDa band seen in Daudi lysate. Recommended for use at 1-3µg/ml.