

GOAT ANTI-RIPK3 / RIP3 ANTIBODY

SKU: EB05004



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	RIP-like protein kinase 3 receptor-interacting protein 3 RIP3 receptor-interacting serine-threonine kinase 3 RIPK3
Accession ID	NP_006862.2
Immunogen	Peptide with sequence C-AWSRPQGWYNHSGK, from the C Terminus of the protein sequence according to NP_006862.2.
Peptide Sequence	C-AWSRPQGWYNHSGK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	11035
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: Approx 45-50kDa band seen in HeLa, Jurkat and NSO (mouse) cell lysates. Recommended for use at 1-3 µg/ml. Overnight incubation is recommended with this antibody.
Application Type	Pep-ELISA