

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RIPK3 / RIP3 ANTIBODY

SKU: EB05004

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias RIP-like protein kinase 3|receptor-interacting protein 3|RIP3|receptor-interacting serine-threonine kinase 3|RIPK3

Names

Accession NP 006862.2

Peptide with sequence C-AWSRPQGWYNHSGK, from the C Terminus of the protein sequence according to

Immunogen NP_006862.2.

Peptide C-AWSRPQGWYNHSGK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Refrigerated

Instructions

Predicted Human **Species**

Human 11035 Gene ID

Product $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/aspiring_medium.png$ Grade

Western Western blot: Approx 45-50kDa band seen in HeLa, Jurkat and NSO (mouse) cell lysates. Recommended for use at

1-3 µg/ml. Overnight incubation is recommended with this antibody. Blot

Application Pep-ELISA **Type**