

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RNF13 / RZF ANTIBODY

SKU: EB08552





Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa	Α	В
50kDa 37kDa		
25kDa		
20kDa		
15kDa	14	





Telephone: (650) 697-3600

SPECIFICATIONS

Unit Size 10

100 μg

Synonym /

Alias RNF13|ring finger protein 13|MGC13689|RZF|RING zinc finger protein

Names

Accession

NP_009213.1; NP_899237.1; NP_899239.2; NP_899240.2

Blocking

Peptide EBP08552

Peptide with sequence HKFKKGDEYDVC, from the internal region of the protein sequence according to Immunogen

NP_009213.1; NP_899237.1; NP_899239.2; NP_899240.2.

Product Comments

This antibody is expected to recognize all reported isoforms (NP_009213.1, NP_899239.2; NP_899240.2).

Peptide

HKFKKGDEYDVC

Sequence Shipping

Shipping Instructions Refrigerated

Predicted

Species

Human, Mouse, Rat, Dog, Cow

Reactive

Species Mouse, Rat

Human

Gene ID

Product

Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Antibody detection limit dilution 1:2000.

Limit

Western Approx 45kDa band observed in Mouse Eye and in Rat Retina lysates (calculated MW of 42.8kDa according to

Blot NP_009213.1). Recommended concentration: 2-4μg/ml.

Application

Type Pep-ELISA, WB

GALLERY IMAGES





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

