


GOAT ANTI-RNF13 / RZF ANTIBODY

SKU: EB08552



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	RNF13 ring finger protein 13 MGC13689 RZF RING zinc finger protein
Names	
Accession ID	NP_009213.1; NP_899237.1; NP_899239.2; NP_899240.2
Blocking Peptide	EBP08552
Immunogen	Peptide with sequence HKFKKGDEYDVC, from the internal region of the protein sequence according to NP_009213.1; NP_899237.1; NP_899239.2; NP_899240.2.
Product Comments	This antibody is expected to recognize all reported isoforms (NP_009213.1, NP_899239.2; NP_899240.2).
Peptide Sequence	HKFKKGDEYDVC
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Mouse, Rat
Human Gene ID	11342
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx 45kDa band observed in Mouse Eye and in Rat Retina lysates (calculated MW of 42.8kDa according to NP_009213.1). Recommended concentration: 2-4µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

