Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RNF25 / AO7 ANTIBODY

SKU: EB05875

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias AO7|FLJ13906|ring finger protein 25|RNF25

Names

Accession NP 071898.2 ID

Blocking

EBP05875 **Peptide**

Peptide with sequence C-TPGSSYPRLPRGQ, from the C Terminus of the protein sequence according to Immunogen

NP 071898.2.

Peptide Sequence

C-TPGSSYPRLPRGQ

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human, Rat, Dog

Human

64320 Gene ID

Rat Gene ID 301515

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Western Blot: No signal obtained yet but low background observed in Human Kidney and Human Brain lysates at Western up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with

Blot other antibodies/lysates?

Application

Pep-ELISA **Type**





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600